Product datasheet Anti-BAG2 Antibody (Clone#AEOE-2) Catalog Number: BM5479

BOSTER®

antibody and ELISA experts
BOSTER BIOLOGICAL TECHNOLOGY

Building C21, 3rd to 5th Floors, Optics Valley Biopharmaceutical Accelerator, East Lake High-Tech Development Zone, Wuhan.

Web: www.boster.com Phone: 027-67845390/1/2 Email: boster@boster.com

Product Name	Anti-BAG2 Antibody (Clone#AEOE-2)	
Gene Name	BAG2	
Source	Rabbit	
Clonality	Monoclonal	
Isotype	IgG	
Species Reactivity	human, mouse, rat	
Tested Application	WB, ICC/IF	
Contents	500 ug/ml; Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide, 0.4-0.5 mg/ml BSA and 50% glycerol.	
Immunogen	A synthesized peptide derived from human BAG2 Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release.	
Concentration	500 ug/ml	
Purification	Affinity-chromatography	
Observed MW	24 kDa	
Dilution Ratios	Western blot (WB):	1:500-2000

Immunocytochemistry/Immunofluorescence (ICC/IF):1:50-200

Storage

12 months from date of receipt, -20°C as supplied.

Background Information

BAG family molecular chaperone regulator 2 is a protein that in humans is encoded by the BAG2 gene. The predicted BAG2 protein contains 211 amino acids. The BAG domains of BAG1, BAG2, and BAG3 interact specifically with the Hsc70 ATPase domain in vitro and in mammalian cells. All 3 proteins bind with high affinity to the ATPase domain of Hsc70 and inhibit its chaperone activity in a Hip-repressible manner. The functional antagonisms displayed between BAG family proteins and Hip suggest that a proper balance of these 2 types of protein is required for achieving optimal cycles of substrate binding and release required for inducting conformational changes in proteins, with Hip promoting peptide substrate binding by Hsc70/Hsp70 and BAG family proteins promoting dissociation.

Selected Validation Data

Product datasheet

Anti-BAG2 Antibody (Clone#AEOE-2)

Catalog Number: BM5479



Building C21, 3rd to 5th Floors, Optics Valley Biopharmaceutical Accelerator, East Lake High-Tech Development Zone, Wuhan.

Web: www.boster.com Phone: 027-67845390/1/2 Email: boster@boster.com

