Product datasheet
Anti-ZO-2/TJP2 Antibody
Catalog Number: PA1971



Building C21, 3rd to 5th Floors, Optics Valley Biopharmaceutical Accelerator, East Lake High-Tech Development Zone, Wuhan.

Web: www.boster.com Phone: 027-67845390/1/2 Email: boster@boster.com

Basic Inform	iation	
Product Name	Anti-ZO-2/TJP2 Antibody	
Gene Name	TJP2	
Source	Rabbit	
Clonality	Polyclonal	
Isotype	lgG	
Species Reactivity	human, mouse, rat	
Tested Application	WB, IHC, ICC/IF	
Contents	500 ug/ml antibody with PBS, 0.02% NaN3, 1 mg/ml BSA and 50% glycerol.	
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of mouse TJP2, identical to the related rat sequence.	
Concentration	500 ug/ml	
Purification	Immunogen affinity purified.	
Observed MW	150 kDa	
Dilution Ratios	Western blot (WB): Immunohistochemistry (IHC): Immunocytochemistry/Immunofluorescence (ICC/IF): (Boiling the paraffin sections in 10mM citrate buffer,pH6.0,o mins is required for the staining of formalin/paraffin sections determined by end user.	

Storage

12 months from date of receipt, -20°C as supplied.

Background Information

TJP2(Tight Junction Protein 2), also known as Zona Occludens 2 or ZO2, is a protein that in humans is encoded by the TJP2 gene. Tight junction proteins(TJPs) belong to a family of membrane-associated guanylate kinase(MAGUK) homologs that are involved in the organization of epithelial and endothelial intercellular junctions. Duclos et al.(1994) mapped the TJP2 gene telomeric to the Friedreich ataxia critical region on chromosome 9q13-q21. TJP2 lies about 70 kb centromeric to the X123 gene and is transcribed in the centromere-to-telomere direction. Using in vitro assays and immunoprecipitation studies, Itoh et al.(1999) showed that the mouse Tjp1, Tjp2, and Tjp3 PDZ1 domains interacted with the C-terminal cytoplasmic domains of Cldn1 through Cldn8. In the mouse inner ear, Walsh et al.(2010) found that Tjp2 expression decreased rapidly between E16.5 and age 1 week to a level in adult mice that was approximately 50% of the level at birth(P0).


Building C21, 3rd to 5th Floors, Optics Valley Biopharmaceutical Accelerator, East Lake High-Tech Development Zone, Wuhan.

Web: www.boster.com Phone: 027-67845390/1/2 Email: boster@boster.com

Reference

Anti-ZO-2/TJP2 Antibody被引用在1文献中。

Selected Validation Data



Western blot analysis of ZO-2/TJP2 using anti-ZO-2/TJP2 antibody (PA1971). The sample well of each lane was loaded with 30 ug of sample under reducing conditions.

Lane 1: human RT4 whole cell lysates,

Lane 2: human A431 whole cell lysates,

Lane 3: human HepG2 whole cell lysates,

Lane 4: human MCF-7 whole cell lysates,

Lane 5: human U2OS whole cell lysates,


Lane 6: human Hela whole cell lysates,

Lane 7: human PC-3 whole cell lysates,

Lane 8: rat RH-35 whole cell lysates,

Lane 9: mouse Hepa1/6 whole cell lysates.

After electrophoresis, proteins were transferred to a membrane. Then the membrane was incubated with rabbit anti-ZO-2/TJP2 antigen affinity purified polyclonal antibody (PA1971) at a dilution of 1:1000 and probed with a goat anti-rabbit IgG-HRP secondary antibody (Catalog # BA1054). The signal is developed using ECL Plus Western Blotting Substrate (Catalog # AR1197). A specific band was detected for ZO-2/TJP2 at approximately 150 kDa. The expected band size for ZO-2/TJP2 is at 131 kDa.



IHC analysis of ZO-2/TJP2 using anti-ZO-2/TJP2 antibody (PA1971) . ZO-2/TJP2 was detected in a paraffin-embedded section of mouse liver tissue. The tissue section was incubated with rabbit anti-ZO-2/TJP2 Antibody (PA1971) at a dilution of 1:200 and developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB (Catalog # AR1027) as the chromogen.


Product datasheet Anti-ZO-2/TJP2 Antibody Catalog Number: PA1971



BOSTER BIOLOGICAL TECHNOLOGY

Building C21, 3rd to 5th Floors, Optics Valley Biopharmaceutical Accelerator, East Lake High-Tech Development Zone, Wuhan.

Web: www.boster.com Phone: 027-67845390/1/2 Email: boster@boster.com



IF analysis of ZO-2/TJP2 using anti-ZO-2/TJP2 antibody (PA1971). ZO-2/TJP2 was detected in an immunocytochemical section of U2OS cells. The section was incubated with rabbit anti-ZO-2/TJP2 Antibody (PA1971) at a dilution of 1:100. DyLight®488 Conjugated Goat Anti-Rabbit IgG (Green) (Catalog # BA1127) was used as secondary antibody. The section was counterstained with DAPI (Catalog # AR1176) (Blue).